OLYMPUS[®]


MagnaFire and MagnaFire SP®

Biomedical Research Cameras

Olympus is about life. About photographic innovations that capture precious moments of life. About advanced medical technology that saves lives. About information and industry-related products that make possible a better living. About adding to the richness and quality of life for everyone. Olympus. Quality products with a **FOCUS ON LIFE**

High Resolution Digital Microscopy Made Easy








Features:

- Full resolution, real time focusing and framing via a Firewire™ interface.
- Live mode features 10-24 images per second for rapid refresh rate and ease of focus.
- Focusing window includes **Auto White Balance**, **Auto Exposure** and Camera Controls.
- Live window features precise frame timing ensuring a **flicker free image** even at short exposures.
- **6-position filter wheel** with red, green, blue, IR-cut and open positions permit a wide spectrum of imaging from UV to IR.
- **Dichromatic, non-polarizing color filters** allows MagnaFire's sensitivity to be at least 12 times greater than cameras using liquid crystal color filter technology.
- FirewireTM interface permits use with laptops (PCMCIA IEEE 1394 OHCI compliant card required).
- Turbo[™] Gain Control reduces exposure duration; eliminates the requirement for binning while maintaining full image resolution at a wide range of exposures.
- Real time image processing within the live window includes: brightness, intensity detect, contrast, histogram equalize, stretch intensity, manual exposure and Turbo™ gain control.
- Color Merge permits individual, user-selectable, image file overlays for imaging with external cubes, filters, and custom mono-, bi-, and tri-chromatic imaging.
- Peltier cooling system combined with a passive radiator maintains the CCD at -40°C ambient, reducing dark current noise and increasing sensitivity for low light level imaging.

Included Software:

- LucisTM DHP; unique contrast enhancement software.
- Image Pro Express™; provides full camera control, image acquisition and processing with database archiving.
- MagnaFire®; provides camera control and image acquisition.



Magnafire SP®

Features:

- **Full resolution,** real time focusing and framing via a Firewire[™] interface.
- Full color, high resolution images with a **single exposure**.
- Automatic "Live" mode, via Firewire™, with Auto White Balance, Auto Exposure. (Speed depending on processor used)
- Live window features precise frame timing ensuring a **flicker free image** even at short exposures.
- Single CCD matrix color sensor produces full color live preview.
- Firewire[™] interface permits use with laptops (PCMCIA IEEE 1394 OHCI compliant card required).
- Turbo[™] Gain Control reduces exposure duration; eliminates the requirement for binning while maintaining full image resolution at a wide range of exposures.
- **Real time image processing** includes: brightness, intensity detect, contrast, histogram equalize, stretch intensity, manual exposure, Turbo[™] gain control.
- Color Merge permits individual, user-selectable, image file overlays for imaging with external cubes, filters, and custom mono-, bi-, or tri-chromatic imaging.
- **Exposure times** from 0.00017 sec to over 10 minutes

Included Software:

- Lucis™ DHP; unique contrast enhancement software.
- Image Pro Express™; provides full camera control, image acquisition and processing with database archiving.
- MagnaFire®; provides camera control and image acquisition.

MagnaFire™ Specifications

MagnaFire Camera Description Megapixel imaging camera system designed for scientific imaging applications

Compatible with Olympus MicroSuite Software

Camera Physical Dimensions 166mm high, 166mm wide, 76mm deep;

weight: 1.6kg

Camera Operating Environment Operating temperature +5°C to +35°C

Storage temperature -10°C to +60°C Humidity 0 to 80% relative non-condensing Shock and vibration: as per ISTA guidelines

Camera Power Supply Power supply dimensions: 213mm long, 155mm

> wide, 64mm tall; weight: 1.7kg Line voltage: 100 - 240V~, 50/60Hz, 33VA;

250mA max current

Fuse type: 5 x 20mm, 400mA, SB, IEC 127

Regulatory Compliance UL/CSA, CE Laboratory Rated

Image Sensor Sony ICX085AL Interline Transfer CCD, 2/3",

1300 x 1030 x 6.7μ square pixels

1280 x 1024 pixels Image Size

Dynamic Range 60dB Signal to Noise Ratio

Color Image Formation Multi-exposure via 6-position filter wheel with

multiple dichroic optical glass color filters. Filters include red, green, blue, IR cut,

clear, black

Filter Optical Performance Transmission: greater than 80% at peak

File Bit Depth 24, 32, 48 bit color: 8 or 16 bit monochrome

Live Image Display Rate 10 images/sec monochrome; 1280 x 1024

IEEE 1394 FireWire™, 400 megabits/sec Computer Interface

 $(Apple^{{\scriptscriptstyle TM}}\ Computer)$

Bright Field Image Capture Time Full field, full resolution color:

less than 4.0 sec.

Full field, full resolution monochrome:

less than 5.0 sec.

CCD Cooling Sealed gas, Peltier cooled to 40°C below

ambient; no desiccant

Chassis Cooling Passive fin radiator

Dark Current 4 electrons/pixel/second @ 0°C nominal

Read Noise 16 electrons/pixel Well Depth 16,000 electrons

Exposure Control Full automatic, manual override

0.00017 sec to 1560.0 sec/color in normal mode. 0.00017~sec to 120.0~sec/color in $Turbo^{\mbox{\tiny TM}}$ mode Manual settings at 1/16 stop intervals

Electronic, software controlled

Camera Shutter

Optical Interface Industry standard C-mount (17.5mm backfocus). No relay lens optics required

Camera Control Windows™ 98 (Second Edition), Windows™ and Display Software 2000. Macintosh™ OS 8.6 or later Microsoft[™] Corporation, Apple[™] Computer

Included Software Lucis[™] DHP; unique contrast enhancement

software

 Image Pro Express™; provides full camera control, image acquisition and processing with database archiving

• MagnaFire®; provides camera control and image acquisition

MagnaFire SP™ Specifications

MagnaFire Camera Description Megapixel color imaging camera system

designed for scientific microscope imaging

applications

Compatible with Olympus MicroSuite Software

Camera Physical Dimensions 47mm high, 166mm wide, 166mm deep;

weight: 1.5kg

Camera Operating Environment Operating temperature +5°C to +35°C

Storage temperature -10°C to +60°C Humidity 0 to 80% relative non-condensing Shock and vibration: as per ISTA guidelines

Camera Power Supply Power supply dimensions: 213mm long, 155mm

wide, 64mm tall; weight: 1.7kg

Line voltage: 100 - 240V~, 50/60Hz, 33VA;

250mA max current

Fuse type: 5 x 20mm, 400mA, SB, IEC 127

Regulatory Compliance UL/CSA, CE Laboratory Rated

Sony ICX085AK Color Interline Transfer CCD, Image Sensor

2/3", 1300 x 1030 x 6.7µ square pixels

Image Size 1280 x 1024 pixels

Dynamic Range 60dB Signal to Noise Ratio 60dB

Color Image Formation Single pass matrix color filter

File Bit Depth 24, 32, 48 bit color: 8 or 16 bit monochrome

Live Image Display Rate 10 images/sec color; 1280 x 1024 resolution,

Computer Interface IEEE 1394 FireWire™, 400 megabits/sec

 $(Apple^{TM} Computer)$

Bright Field Image Capture Time Full field, full resolution color:

less than 2.0 sec.

Read Noise 16 electrons/pixel

Well Depth 16,000 electrons

Exposure Control Full automatic, manual override

0.00017 sec to 10+ min./color in normal mode 0.00017~sec to 1~min./color in $Turbo^{\mbox{\tiny TM}}$ mode Manual settings at 1/16 stop intervals

Camera Shutter Electronic, software controlled

Optical Interface Industry standard C-mount (17.5mm backfocus). No relay lens optics required

Windows™ 98 (Second Edition), Camera Control

and Display Software Windows™ 2000


Included Software \bullet Lucis $^{\scriptscriptstyle{TM}}$ DHP; unique contrast enhancement

software

 Image Pro Express™; provides full camera control, image acquisition and processing with

database archiving

· MagnaFire®; provides camera control and image acquisition



DLYMPUS

Focus On Life

1-800-446-5967 • www.olympus.com

Olympus America Inc. Scientific Equipment Group Two Corporate Center Drive Melville, NY 11747-3157 800-446-5967 Fax: 631-844-5112

Canada: Carsen Group Inc. 151 Telson Road Markham, Ontario L3R 1E7 905-479-4100 Fax: 905-479-2595

Latin America: Olympus America Inc. Latin America Group 6100 Blue Lagoon Drive, Suite 390 Miami, FL 33126-2067 305-266-2332 Fax: 305-266-9168

Specifications subject to change without notice.